



# ePrep Location

Pub No. 99-10010-01




Place ePrep on a clean flat surface with no obstruction to the axis movement. Adjust ePrep feet to level instrument if required.

## Connect USB Cable and Open Tablet Software



Power up the Tablet Computer and open the ePrep software via the desktop ePrepApp icon shown above. Further installation instructions can be found on the Tablet.

# ePrep Software Main Menu



Go to System>Quick Start for further setup instructions.

# IMPORTANT SAFETY WARNINGS

---



**Do not enter the instrument work space or place any part of your body near moving axes during operation.**

---



If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired. Good Laboratory Practice (GLP) should be followed when using the ePrep. Refer: 21CFR58, 21CFR58, Directive 2004/9/EC and specific product GLP Directives

---



ePrep is fitted with an Emergency Stop Button found at the front right side of the instrument. Pressing the button will cut off all power to the instrument. To reinstate power, rotate the emergency button 90° clockwise.


---



Only approved accessories and consumables can be used with the instrument. Eprep assumes no liability for the customer's failure to comply with these requirements.

---

# Install Drip Tray



Drip Tray Position

Install the Drip Tray below the deck area. The Tray is located in the front panel (first removed panel) of the ePrep box.

## Manually Park the System




Manually move the  
XCHANGE axis onto  
the Park Station

Park Station


Before powering up the ePrep, manually move axis so XCHANGE change coupling is sitting on the Park Station.

## Wash Station Electrical Connection



Ensure ePrep power is switched off. Before installing the Wash Station deck plate, connect Pumps to ePrep via the electrical connection shown.


## Install the Wash Station Deck Plate



Install the Wash Station deck plate and lock it down using in installation screw.




## Wash Station Liquid Connections




Install tubing to the Wash (IN) and Wash (OUT). Connect Wash (IN) to a receptacle bottle with at least 600mL of wash solution and wash (OUT) to a waste receptacle bottle.


## Install Rear-Right Plate




Locate and lock the Rear Right deck plate to the ePrep frame using the supplied installation screws.

## Cable Connection



Connect the Power Pack to ePrep and the mains power using a Country Specific IEC C13 cable (NOT SUPPLIED). Connect Tablet to ePrep with the supplied USB cable. Switch ePrep on. When the Tablet has established USB connection with ePrep the  icon will appear at the top right of the tablet screen.

# Emergency Stop Button



Emergency Stop Button


Ensure the Emergency Stop button is released by turning it 90° in the direction of the arrow. Note: press to engage (switches power OFF to ePrep) and rotate to release allows ePrep to be powered by the Power Switch.

# Status Indicator Lights



<b>NO LIGHT</b>	<b>No power to the instrument</b>	<b>ORANGE</b>	Workflow paused by user with the option to abort or continue.
<b>SOLID BLUE</b>	ePrep ready for operation	<b>RED</b>	Bumper or instrument error. Check for obstructions to the rotating foot. Typically, items placed on the deck or different height vial in a rack. Faults with instrument hardware, see trouble shooting information.
<b>GREEN</b>	Workflow running and will remain green until completion or interruption of workflow		
<b>FLASHBLUE</b>	Workflow completed successfully		
<b>FLASH PURPLE</b>	Workflow aborted by user	<b>WHITE</b>	Mains power disconnected and ePrep running on battery power. Battery life approximately 20 min.
<b>FLASH RED</b>	Workflow aborted by system or Failure to detect racks or Failure to detect tubes.		

## Install Vial Adapter Plates for Calibration



Adapter Plates  
in Position 1 and 15


Install supplied Adapter Plates in deck positions 1 and 15. Remove all other Adapter Plates from the deck.

## Calibration Probe



Locate the Calibration Probe found in the ePrep Accessories Box.


## Remove the Sensor Retaining Screw



Use the Hex end of the tool to undo the Sensor Retaining Screw at the base of the ePrep Rotating Foot. Hold the Bumper so that it does not fall off when the screw is removed.



## Install Calibration Probe



Screw end of the  
Calibration Probe

Replace the screw with the Calibration Probe and firmly screw it into the threaded hole, ensuring the foot sensor remains in position. **WARNING:** Do not overtighten.

# Check Installation of Calibration Plates

**ePREP** Calibrate Deck

Back Info User Comm Battery

**Calibration Options**

- Axis calibration
  - X and Y Axis
  - Z Axis
- Tool Rack Calibration
  - Tool Y Position
  - Tool Z position
- Foot Calibration
  - Foot Bumper
  - Foot Rotation

**Calibration Data**

ePrep Serial No	
4	
X Offset	X Grad
2083	0.00625
Y Offset	Y Grad
2436	-0.0006779
Z Offset	Z Grad/Inch
6800	-7.256-05
Z Grad/Inch	Z BackToFromOff
0.0005318	-0.5720215
X Coeff	Y Coeff
0.0014946	-0.0007916
ToolRackOffset Y	ToolRackOffset Z
-30.71701	32.02471
Bumper Offset	Rotate Offset
-1.495319	14
Spare 1	Calib Date
15	161220

Install the calibration Mats at position 1 and 15 and remove all other mats

ARE THE CALIBRATION MATS INSTALLED AND POSITIONED CORRECTLY  
Press Help for installation instructions

HELP Yes No

Calibrate Deck  
Test Calib Points  
Initialise

Check the Adapter Plates are installed correctly. On the software Main Screen click on System > Service > Calibrate to open the Calibrate Deck option. Software will ask if Plates are positioned correctly.

# Check Installation of Calibration Probe

The screenshot shows the 'Calibrate Deck' interface of the ePREP system. The interface includes a top navigation bar with icons for Back, Info, User, Comm, and Battery. On the left, there is a sidebar with 'Calibration Options' (Axis calibration, X and Y Axis, Z Axis), 'Tool Rack Calibration' (Tool Y Position, Tool Z position), 'Foot Calibration' (Foot Bumper, Foot Rotation), and buttons for 'Calibrate Deck', 'Test Calib Points', and 'Initialise'. The main area features a large dark box with the text 'Insert picture of installed probe here' and a confirmation dialog asking 'IS CALIBRATION PROBE INSTALLED CORRECTLY?' with 'Help', 'Yes', and 'No' buttons. A warning message at the bottom states: 'Installation of calibration probe is not required when calibrating: Z axis, Tool-Y, Tool-Z, or Foot Bumper and must be removed when calibrating Foot Rotation'. On the right, the 'Calibration Data' table is displayed.

Calibration Data	
ePrep Serial No 4	
X Offset	X Grad
2083	0.00625
Y Offset	Y Grad
2436	-0.0006779
Z Offset	Z Grad/Rev
6800	-7.25E-05
Z Grad/Rev	Z BackToFrontOff
0.0005318	-0.5720215
X Coeff	Y Coeff
0.0014946	-0.0007916
ToolRackOffset Y	ToolRackOffset Z
-30.71701	32.02471
Bumper Offset	Rotate Offset
-1.495319	14
Spare 1	Calib Date
15	161220

If required the software will ask if the Calibration Probe has been correctly installed.


# Deck Calibration

If Adapter Plates and Probe are installed, the following calibration options are available: X-Axis, Y-Axis, Z-Axis, Tool-Y, Tool-Z, Foot Sensor, Foot Rotation Note: user will be prompted to remove the calibration probe for foot rotation \*.

**Calibration Data** panel records and displays the calibration data used to calculate operational parameters. Calibration points are saved on the instrument firmware for real time calculation of positional data.

**Calibrate Deck:** starts or stops (Aborts) the calibration. If required, a system initialisation will be performed.

**Test Calibrate Points** button can be used to test the X, Y and Z calibration at four points on the deck. If required user will be prompted to remove the probe during the procedure.



# Select ALL Calibration Options


Calibration Options

The screenshot shows the 'Calibrate Deck' interface of the ePREP software. On the left, a sidebar lists calibration options under three categories: 'Axis calibration' (X Axis, Y Axis, Z Axis), 'Tool Rack Calibration' (Tool Y Position, Tool Z position), and 'Foot Calibration' (Foot Bumper, Foot Rotation \*). A yellow bracket groups these options under the label 'Calibration Options'. Below the options are three buttons: 'Calibrate Deck', 'Test Calib Points', and 'Initialise'. The main area features a top video feed showing a deck with four points and a bottom section titled 'Calibration Data' which is currently empty. On the right, a 'Calibration Data' table displays various parameters and their values.

Calibration Data	
ePrep Serial No 4	
X Offset	X Grad 2021 0.00255
Y Offset	Y Grad 2383 -0.0006581
Z Offset	Z GradRear 7122 -0.000161
Z GradFront	Z BackToFrontOff
0.0002863	-0.3230896
X Coeff	Y Coeff
0.0004425	0.0015543
ToolRackOffset Y	ToolRackOffset Z
-27.47032	33.58932
Bumper Offset	Rotate Offset
-0.0465424	-30
Spare 1	Calib Date
15	170110

On initial installation and commissioning, all calibration options should be selected. Note: the software does allow for selective “optional” calibration depending on installation of Probe and Adapter Plates.


# Trigger Calibrate Deck



Calibrate Deck Button

Click on “Calibrate Deck”. ePrep will move the Calibration Probe to allocated positions on the deck and tool rack to record instrument calibration data. Full calibration takes around 10 minutes.

# Calibration Completion and Data Write



When complete the calibration dialogue will display "Calibration updated successfully". Calibration data is automatically written to the instrument EEPROM to be used during operation.

## Remove Calibration Probe




Remove the Probe and  
replace the Holding Screw

When prompted, remove the Calibration Probe and replace the Holding Screw before proceeding. Store the Calibration Probe in the threaded hole next to the Wash Station.




## Store Calibration Probe




Screw the Calibration Probe into its deck position next to the Wash Station for storage.

## Installing Syringes and Tools




ePrep Syringes and tools are fitted with RFID identification collars. Place the tool in the cassettes in the correct orientation (as shown) with the plunger fully down and the Tool Retainer Clip covering the plunger button.

## Vial Adapter Plates and Vials



Vial Adapter Plates are polarised with the rear pins larger than the front pins. The engraved Plate Position Indicator is used to align the rack position number on the front of the deck.

# Software Specified Vial Racks and Tools - IMPORTANT



ONLY Vial Racks, Syringes and Tools that are **specified** in the software can be used with ePrep. On initial identification, Tools must be **registered** to enable the software to identify registered assets have corresponding software profile for simple programming within a Workflow.

## RFID Rack Identification

RFID Tag on Post



RFID Tag on Rack

ePrep requires an RFID tag for identification of all vial racks on the deck. RFID tags are located on the rack itself or on a plate post. All vials within a rack must be the same height and type as the first vial because only the height of the first vial is sensed.